

Recombinant Human Ephrin B Receptor 1,EphB1 (ICD, C-6His)

Catalog No: #AP75388

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75388-1 10ug #AP75388-2 50ug #AP75388-3 500ug #AP75388-4 1mg

Description

Product Name	Recombinant Human Ephrin B Receptor 1,EphB1 (ICD, C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P54762
Uniprot	P54762
GeneID	2047;
Calculated MW	48.8kDa
Target Sequence	MSRKRAYSKEAVYSDKLQHYSTGRGSPGMKIYIDPFTYEDPNEAVREFAKEIDVSFVKIEEVIGAGEFGGEVYKG RLKLP GKREIYVAIKTLKAGYSEKQRRDFLSEASIMGQFDHPNIIRLEGVVTKSRPVMIIITEFMENGALDSFLRQ NDGQFTVIQLVGMLRGIAAGMKYLAEMNYVHRDLAARNILLNSNLVCKVSDFGLSRYLQDDTSDPTYTSSLGG KIPVRWTAPEAIAIYRKFTSASDVWSYGIVMWEVMSFGERPYWDMSNQDVINAIEQDYRLPPPMDCPAALHQL MLDCWQKDRNSRPRFAEIVNTLDKMIRNPASLKT VATITAVPSQPLLD RSIPDFTAFTTVDDWLSAIKMVQYRD SFLTAGFTSLQLVTQMTSEDLLRIGITLAGHQKILNSIHSMRVQISQSPTAMAVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM Tris, 150mM NaCl, pH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only