

Recombinant Human Fumarylacetoacetase,FAH (C-6His)



Catalog No: #AP75389

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75389-1 10ug #AP75389-2 50ug #AP75389-3 500ug #AP75389-4 1mg

Description

Product Name	Recombinant Human Fumarylacetoacetase,FAH (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16930
Uniprot	P16930
GeneID	2184;
Calculated MW	47.4kDa
Target Sequence	SFIPVAEDSDFPIHNLPYGVFSTRGDPRPRIGVAIGDQILDLSIIKHLFTGPVLSKHQDVFNQPTLNSFMGLGQA AWKEARVFLQNLLSVSQARLRDDTELRKCAFISQASATMHLPATIGDYTDFYSSRQHATNVGIMFRDKENALM PNWLHLPVGYHGRASSVVVSGTPIRRPMGQMKPDDSKPPVYGACKLLDMELEMAFFVGPNGRLGEPISKKA HEHIFGMVLMNDWSARDIQWEYVPLGPFGLKSFGTTVSPWVVPMDALMPFAVNPKNQDPRPLPYLCHDEP YTFDINLSVNLKGEQMSQAATICKSNFKYMYWTMLQQLTHHSVNGCNLRPGDLLASGTISGPEPENFGSMLEL SWKGTKPIDLNGQTRKFLLDGDEVIITGYCQGDGYRIGFGQCAGKVL PALLPSVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only