

Recombinant Human Coagulation Factor XIII A Chain (C-6His)



Catalog No: #AP75391

Orders: order@signalwayantibody.com

Package Size: #AP75391-1 10ug #AP75391-2 50ug #AP75391-3 500ug #AP75391-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Coagulation Factor XIII A Chain (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P00488
Uniprot	P00488
GeneID	2162;
Calculated MW	80.3kDa
Target Sequence	GVNLQEFLNVTSVHLFKERWDTNKVDHHTDKYENNKLIVRRGQSFYVQIDFSRPYDPRRD LFRVEYVIGRYPQ ENKGTYPVPIVSELQSGKWGAKIVMREDRSVRLSIQSSPKCIVGKFRMYVAVWTPYGLVLR TSNRPETD TYILF NPWCEDDAVYLDNEKEREEYVLNDIGVIFYGEVNDIKTRSWSYQQFEDGILDTCLYMDRAQMDLSGRGNPIK VSRVGSAMVNAKDDEGLVGSWDNIYAGVPPSAWTGSVDILLEYSSENVPVRYGQCWVFAGVFNTFLRCL GIPARIVTNYFSAHDNDANLQMDIFLEEDGNVNSKLTKDSVWNYHCWNEAWMTRPDLVPGFGGWQAVDSTP QENS DGM YRCGPASVQA IKHGHVCFQFDAPFVFAEVNSDLIYITAKKDGTHVVENVDATHIGKLIVTKQIGGDG MMDITD TYK FQEGQEEERLAL ETALMYGAKKPLNTEGVMKSRSNVDMDFEVENAVLGKDFKLSITFRNNSHN RYTITAYLSANITFYTGVPKAEFKKETFDTLEPLSFKKEAVLIQAGEYMGQLLEQASLHFFVTARINETRDVLAK QKSTVLTIP EIIIKV RGTQVVGSDMTVTVQFTNPLKETLRNVVWVHLDGPGVTRPMKKMFREIRPNSTVQWEEV CRPWVSGHRKLIASMSSDSL RHVYGELDVQIQRPSMVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 50 mM NaCl, 5% Sucrose, 1% Tween 20 (v,v), 0.3% Histidine (w,v), pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only