

Recombinant Human Sorbitol Dehydrogenase, SORD (C-6His)



Catalog No: #AP75393

Orders: order@signalwayantibody.com

Package Size: #AP75393-1 10ug #AP75393-2 50ug #AP75393-3 500ug #AP75393-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Sorbitol Dehydrogenase, SORD (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q00796
Uniprot	Q00796
GeneID	6652;
Calculated MW	39.3kDa
Target Sequence	<p>AAAAPNNLSLVVHGPGDLRLNYPEPGPNEVLLRMHSVGICGSDVHYWEYGRIGNFIVKKPMVLGHEASG TVEKVGSSVKHLKPGDRVAIEPGAPRENDEFCKMGRYNLSPSIFFCATPPDDGNLCRFYKHNAAFCYKLPDN VTFEEGALIEPLSVGIHACRRGGVTLGHKVLVCGAGPIGMVTLVAKAMGAAQVVVTDLSATRLSKAKEIGADL VLQISKESPQEIARKVEGQLGCKPEVTIECTGAEASIQAGIYATRSGGTLVLVGLGSEMTPVLLHAAIREVDIKG VFRYCNTWPVAISMLASKSVNVKPLVTHRFPLEKALEAFETFKKGLGLKIMLKCDPSDQNPVDHHHHHH</p>
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 0.2M NaCl, 5mM DTT, 20% glycerol, pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only