

Recombinant Human Acyl-Coenzyme A Thioesterase 13, ACOT13 (C-6His)



Catalog No: #AP75394

Orders: order@signalwayantibody.com

Package Size: #AP75394-1 10ug #AP75394-2 50ug #AP75394-3 500ug #AP75394-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Acyl-Coenzyme A Thioesterase 13, ACOT13 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NPJ3
Uniprot	Q9NPJ3
GeneID	55856;
Calculated MW	15.9kDa
Target Sequence	TSMTQSLREVIKAMTKARNFERVLGKITLVSAAPGKVICEMKVEEHTNAIGTLHGGLTATLVDNISTMALLCTE RGAPGVSVDNMNITYMSPAKLGEDIVITAHVLKQGKTLAFTSVDLTNKATGKLIAQGRHTKHLGNVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only