

Recombinant Human GALNTL1 (C-6His)

Catalog No: #AP75399



Package Size: #AP75399-1 10ug #AP75399-2 50ug #AP75399-3 500ug #AP75399-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GALNTL1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8N428
Uniprot	Q8N428
GeneID	57452;
Calculated MW	61kDa
Target Sequence	DNRAHAASSGGRGAQRAGRRSEQLREDRTIPLIVTGTPSKGFDEKAYLSAKQLKAGEDPYRQHAFNQLESDK LSPDRPIRDTRHYSCPSVSYSSDLPATSVIITFHNEARSTLLRRTVKSVLNRTAPANLIQEIILVDDFSSDPEDCLLLT RIPKVKCLRNDRREGLIRSRVRGADVAAATVLTFLDSHCEVNTWELPPMLQRVKEDHTRVVSPIIDVISLDNFA YLAASADLRGGFDWSLHFKWEQIPLEQKMTRTDPTRPIRTPVIAGGIFVIDKSWFNHLGKYDAQMDIWGGENF ELSFVRVVMCGGSLEIVPCSRVGHVFRKRHPYNFPEGNALT YIRNTKRTAEVWMDEYKQYYEARPSAIGKAF GSVATRIEQRKKNCKSFRWYLENVPELTVPVKEALPGIIKQGVNCLSQGQNTAGDFLLGMGICRGSKAKNP QPAQAWLFSDHLIQQQGKCLAATSTLMSSPGSPVILQMCNPREGKQKWRKGSFIQHSVSGLCLETKPAQLV TSKCQADAQAQQWQLLPHTVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only