

Recombinant Human Membrane Primary Amine Oxidase,AOC3(C-6His)

Catalog No: #AP75407

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75407-1 10ug #AP75407-2 50ug #AP75407-3 500ug #AP75407-4 1mg

Description

Product Name	Recombinant Human Membrane Primary Amine Oxidase,AOC3(C-6His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q16853
Uniprot	Q16853
GeneID	8639;
Calculated MW	82.6kDa
Target Sequence	<p> RGGDGGEPSQLPHCPSPVSPSAQPWTHPGQSQLFADLSREELTAVMRFLTQRLGPGLVDAQAARPSDNCVF SVELQLPPKAAALAHLDRGSPPPAREALAIFFGRQPQPNVSELVVGPLPHPSYMRDVTVERHGGPLPYHRR PVLFAQEYLDIDQMIFNRELPQASGLLHCCFYKHRGRNLVTMTTAPRGLQSGDRATWFGLYYNISGAGFFLHH VGLELLVNHKALDPAWRTIQKVFYQGRYYDSLAQLEAQFEAGLVNVVLIIPDNGTGGSSWSLKSPVPPGPAPPL QFYPPQGRFVSQGSRVASSLWTFSGFLGAFSGPRIFDVRFQGERLVEISLQEALAIYGGNSPAAMTTRYVD GGFGMGKYTTPLTRGVDCPYLATYVDWHFLLESQAPKTIRDAFCVFEQNQGLPLRRHSDLYSHYFGGLAET VLVVRSMTLLNYDYVWDTVFHPSGAIEIRFYATGYISSAFLFGATGKYGNQVSEHTLGTVHTSAHFKVDLDV AGLNWWAEDMVFPMAVPWSPEHQLQRLQVTRKLEEMEEQAFLVGSATPRYLYLASNHSNKWGHPRG YRIQMLSFAGEPLPQNSSMARGFSWERYQLAVTQRKEEEPSSSVFNQNDPWAPTVDVDFSDFINNETIAGKDL VAWVTAGFLHIPHAEDIPNTVTVGNGVGFLLRPYNFFDEDEPSFYSADSIYFRGDQDAGACEVNPLACLQAAA CAPDLPAFSGGGFSHNHHHHHH </p>
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris, 500mM NaCl, pH8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only