

## Recombinant Human Serpin A4, Kallistatin (C-6His)

Catalog No: #AP75408



Package Size: #AP75408-1 10ug #AP75408-2 50ug #AP75408-3 500ug #AP75408-4 1mg

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                 |  |
|-----------------|--|
| Product Name    | Recombinant Human Serpin A4, Kallistatin (C-6His)  |
| Host Species    | Human Cells  |
| Purification    | Greater than 95% as determined by reducing SDS-PAGE.   |
| Specificity     | Human  |
| Accession No.   | P29622   |
| Uniprot         | P29622   |
| GeneID          | 5267;  |
| Calculated MW   | 47.1kDa  |
| Target Sequence | QLHVEHDGESCNSSSHQILETGEGSPSLKIAPANADFAFRFYLLIASETPGKNIFFSPLSISAAYAMLSLGACSHRSQILEGLGFNLTESESDVHRGFQHLLHTLNLPGHGLETRVGSALFLSHNLKFLAKFLNDTMAVYEAKLFHTNFYDTVGTIQLINDHVKKETRGKIVDLVSELKDKDVLMLVNYIYFKALWEKPFISSRTTPKDFYVDENTTVRVP<br>MMLQDQEHHWYLDHRYLPCSVLRMDYKGDATVFFILPNQGKMREIEEVLTPEMPLRWNNLLRKRNFYKKLELHLPKFSISGSYVLDQILPRLGFTDLFSKWADLSGITKQKLEASKSFHKATLDVDEAGTEAAAATSFAIKFFSAQ<br>TNRHILRFNRPFLLVIFSTSTQSVLFLGKVVDPKPHHHHHH |
| Formulation     | Lyophilized from a 0.2 µm filtered solution of 50mM TrisHCl, 150mM NaCl, 10mM NaCl, pH 8.0.  |
| Storage         | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.   |

Note: This product is for in vitro research use only