

# Recombinant Human HDGFRP3,HRP-3,HDGF2 (C-6His)

Catalog No: #AP75415

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75415-1 10ug #AP75415-2 50ug #AP75415-3 500ug #AP75415-4 1mg

## Description

Product Name	Recombinant Human HDGFRP3,HRP-3,HDGF2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y3E1
Uniprot	Q9Y3E1
GeneID	50810;
Calculated MW	23.6kDa
Target Sequence	MARPPREYKAGDLVFAKMKGYPHWPARIDELPEGAVKPPANKYPIFFFGTHETAFLGPKDLFPYKEYKDKFG KSNKRKGFNEGLWEIENNPVKFTGYQAIQQSSSETEGEGGNTADASSEEEGDRVEEDGKGRKNEKAGS KRKKSYSKSSKQSRKSPGDEDDKDCKEEENKSSSEGGDAGNDTRNTTSDLQKTSEGTVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM MES, 150mM NaCl, pH5.5.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only