

Recombinant Human Ameloblastin, AMBN (C-6His)

Catalog No: #AP75418



Package Size: #AP75418-1 10ug #AP75418-2 50ug #AP75418-3 500ug #AP75418-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------|---|
| Product Name | Recombinant Human Ameloblastin, AMBN (C-6His) |
| Host Species | Human Cells |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | Q9NP70 |
| Uniprot | Q9NP70 |
| GeneID | 258; |
| Calculated MW | 46.4kDa |
| Target Sequence | VPFFPQQSGTPGMASLSLETMRQLGSLQRLNTLSQYSRYGFGKSFNSLWMHGLLPPHSSLPWMRPREHET QQYEYSLPVHPPPLPSQPSLKPQQPGLKPFLLQSAATTNQATALKEALQPPHILGHLPLQEGELPLVQQQVAP SDKPPKPELPGVDFADPQGPSPLPGMDFDPQGPSPLGLDFADPQGSTIFQIARLISHGMPQNKQSPLYPGM LYVPFGANQLNAPARLGIMSSEEVAGGREDPMAYGAMFPGFGGMRPGFEGMPHNPAMGGDFTLEFDSPVA ATKGPENEEGGAQGSPMEANPDNLENPAFLTELEPAPHAGLPALPKDDIPGLPRSPSGKMKGLPSVTPAAA DPLMTEPELADVYRTYDADMTTSVDFQEEATMDTTMAPNSLQTSMPGNKAQEPMMHDAWHFQEPVDHHHH HH |
| Formulation | Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4. |
| Storage | Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months. |

Note: This product is for in vitro research use only