

Recombinant Human IL-22 Receptor Subunit α2,IL-22BP,IL-22RA2 (C-6His)

Catalog No: #AP75420

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75420-1 10ug #AP75420-2 50ug #AP75420-3 500ug #AP75420-4 1mg

Description

Product Name	Recombinant Human IL-22 Receptor Subunit α2,IL-22BP,IL-22RA2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q969J5-2
Uniprot	Q969J5
GeneID	116379;
Calculated MW	25.7kDa
Target Sequence	TQSTHESLKPQRVQFQSRNFHNLQWQPGRALTGNSSVYFVQYKIYGQRQWKNKEDCWGTQELSCDLTSET SDIQEPYYGRVRAASAGSYSEWSMTPRFTPWWETKIDPPVMNITQVNGSLLVILHAPNLPYRYQKEKNVSIED YYELLYRVFIINNSLEKEQKVYEGAHRAVEIEALTPHSSYCVVAEIQPMLDRRSQRSEERCVEIPVDHHHHHHH
Formulation	Lyophilized from a 0.2 μm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only