

Recombinant Human Pyruvate Kinase, Liver And RBC,PKLR (C-6His)



Catalog No: #AP75422

Orders: order@signalwayantibody.com

Package Size: #AP75422-1 10ug #AP75422-2 50ug #AP75422-3 500ug #AP75422-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Pyruvate Kinase, Liver And RBC,PKLR (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P30613
Uniprot	P30613
GeneID	5313;
Calculated MW	62.9kDa
Target Sequence	MSIQENISSLQLRSWVSKSQRD LAKSILIGAPGGPAGYLRRASVAQLTQELGTAFFQQQLPAAMADTFLEHL CLLDIDSEPVAARSTSIATIGPASRSVERLKEMIKAGMNIARLNFSHGSHEYHAESIANVREAVESFAGSPLSYR PVAIALDTKGPEIRTGILQGGPESEVELVKGSQVLVTVDPFRTRGNANTVWVDYPNIVRVVPVGGRIYIDDGLI SLVVQKIGPEGLVTQVENGGLVLSRKGVNLPGAQVDLPGLSEQDVRDLRFGVEHGVDFVASFVRKASDVAA VRAALGPEGHGIKIISKIENHEGVKRFDEILEVSDGIMVARGDLGIEIPAQVFLAQKMMIGRCNLGKPVVCAT QMLESMITKPRPTRAETSDVANAVLDGADCIMLSGETAKGNFPVEAVKMQHAIAREAEAAVYHRQLFEELRRA APLSRDPTEVTAIGAVEAAFKCCAAAIIVLTTTGRSAQLLSRYRPRAAVIAVTRSAQAARQVHLCRGVFPFLYRE PPEAIWADDVDRRVQFGIESGKLRGFLRVGDLVIVVTGWRPGSGYTNIMRVLISVDHHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM Tris-HCl, 1mM DTT 0.2M NaCl, 10% glycerol, pH 8.0.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only