

Recombinant Human Apolipoprotein A-IV,ApoA4 (C-6His)

Catalog No: #AP75436

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75436-1 10ug #AP75436-2 50ug #AP75436-3 500ug #AP75436-4 1mg

Description

Product Name	Recombinant Human Apolipoprotein A-IV,ApoA4 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P06727
Uniprot	P06727
Calculated MW	44.44kDa
Target Sequence	EVSAQVATVMWDYFSQLSNNAKEAVEHLQKSELTQQLNALFQDKLGEVNTYAGDLQKKLVPFATELHERLA KDSEKLKEEIGKELEELRARLLPHANEVSQKIGDNLRELQQRLEPYADQLRTQVNTQAEQLRRQLTPYAQRME RVLRENADSLQASLRPHADELKAKIDQNVEELKGRLLPYADEFKVKIDQTVELRRSLAPYAQDTQEKLNHQL EGLTFQMKNNAEELKARISASAEELRQLAPLAEDVRGNLRGNTEGLQKSLAELGGHLDQQVEEFRRRVEPY GENFNKALVQQMEQLRQKLGPHAGDVEGHLSEFLKDLRDKVNSFFSTFKKESQDKTSLPELEQQQEQQQ EQQQEQVQMLAPLESVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH 7.2.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only