

Syk (Phospho-Tyr348) Polyclonal Antibody

Catalog No: #12293

Package Size: #12293-1 50ul #12293-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Syk (Phospho-Tyr348) Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Hu Ms Rt Mk
Specificity	Phospho-Syk (Y348) Polyclonal Antibody detects endogenous levels of Syk protein only when phosphorylated at Y348.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Syk around the phosphorylation site of Y348.
Target Name	Syk
Modification	Phospho
Other Names	SYK; Tyrosine-protein kinase SYK; Spleen tyrosine kinase; p72-Syk
Accession No.	Swiss-Prot: P43405NCBI Gene ID: 6850
Uniprot	P43405
GeneID	6850;
Target Species	human
SDS-PAGE MW	72kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1/500 - 1/2000

Immunohistochemistry: 1/100 - 1/300


ELISA: 1/10000

Not yet tested in other applications

Images



Western Blot analysis of 3T3 cells using Phospho-Syk (Y348) Polyclonal Antibody



Western Blot analysis of NIH-3T3 cells using Phospho-Syk (Y348) Polyclonal Antibody

Note: This product is for in vitro research use only