

Recombinant Human Follitropin Subunit β ,FSHB (C-6His)

Catalog No: #AP75441

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75441-1 10ug #AP75441-2 50ug #AP75441-3 500ug #AP75441-4 1mg

Description

Product Name	Recombinant Human Follitropin Subunit β ,FSHB (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P01225
Uniprot	P01225
GeneID	2488;
Calculated MW	13.5kDa
Target Sequence	NSCELTNITIAIEKEEFCISINTTWCAGYCYTRDLVYKDPARPKIQKTCTFKELVYETVRVPGCAHHADSLYT YPVATQCHCGKCDSSTDCTVRGLGPSYCSFGEMKEVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only