

Recombinant Human N-Sulphoglucosamine Sulphohydrolase,SGSH (C-6His)



Catalog No: #AP75442

Orders: order@signalwayantibody.com

Package Size: #AP75442-1 10ug #AP75442-2 50ug #AP75442-3 500ug #AP75442-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human N-Sulphoglucosamine Sulphohydrolase,SGSH (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P51688
Uniprot	P51688
GeneID	6448;
Calculated MW	55.72kDa
Target Sequence	RPRNALLLLADDGGFESGAYNNSAIATPHLDALARRSLLFRNAFTSVSSCSPSRASLLTGLPQHONGMYGLHQ DVHFNFSFDKVRSLPLLLSQAGVRTGIIGKKHVGPEVYPFDFAYTEENGSVLQVGRNITRIKLLVRKFLQTQD DRPFFLYVAFHDPHRCGHSQPQYGTFCFKFNGESGMGRIPDWTPQAYDPLDVLVPYFVPNTPAARADLAA QYTTVGRMDQGVGLVLQELRDAGVLNDTLVIFTSDNGIPFSGRTNLYWPGTAEPLLVSSPEHPKRWGQVSE AYVSLDLTPTILDWFSIPYPSYAIFGSKTIHLTGRSLLPALEAEPLWATVFGSQSHHEVTMSYPMRSVQHRHF RLVHNLNFKMPFPIDQDFYVSPFQDLLNRTTAGQPTGWYKDLRHYYRARWELYDRSRDPHETQNLATDPR FAQLEMLRDQLAKWQWETHDPWVCAPDGVLEEKLSPPQCQPLHNELVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM GaCl ₂ , 10%Glycerol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only