

Recombinant Human Ribonuclease Pancreatic,RNASE1 (C-6His)



Catalog No: #AP75455

Orders: order@signalwayantibody.com

Package Size: #AP75455-1 10ug #AP75455-2 50ug #AP75455-3 500ug #AP75455-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Ribonuclease Pancreatic,RNASE1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P07998
Uniprot	P07998
GeneID	6035;
Calculated MW	15.6kDa
Target Sequence	KESRAKKFQRQHMDSDSSPSSSSTYCNQMMRRRNMTQGRCKPVNTFVHEPLVDVQNVCFQEKVTCKNGQ GNCYKSNSSMHITDCRLTNGSRYPNCAYRTSPKERHIIIVACEGSPYVPVHFDASVEDSTVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, 10% Glycerol, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only