

Recombinant Human Chymotrypsin-Like Elastase Family Member 3A,CELA3A (C-6His)



Catalog No: #AP75456

Orders: order@signalwayantibody.com

Package Size: #AP75456-1 10ug #AP75456-2 50ug #AP75456-3 500ug #AP75456-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Chymotrypsin-Like Elastase Family Member 3A,CELA3A (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09093
Uniprot	P09093
GeneID	10136;
Calculated MW	28.9kDa
Target Sequence	SGYGPPSSHSSSRVVHGEDAVPYSWPWQVSLQYEKSGSFYHTCGGSLIAPDWVVTAGHCISRDLTYQVVLG EYNLAVKEGPEQVIPINSEELFVHPLWNRSCVACGNDIALIKLSRSAQLGDAVQLASLPPAGDILPNKTPCYITG WGRLYTNGPLPDKLQARLPVVDYKHCSRWNWWGSTVKKTMVCAGGYIRSGCNGDSSGGLNCP TEDGGW QVHGVTSFVSGFGCNFIWKPTVFTRVSAFIDWIEETIASHVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,10%Glycerol,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only