

Recombinant Human Carboxypeptidase B2,CPB2 (C-6His)

Catalog No: #AP75457

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75457-1 10ug #AP75457-2 50ug #AP75457-3 500ug #AP75457-4 1mg

Description

Product Name	Recombinant Human Carboxypeptidase B2,CPB2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96IY4
Uniprot	Q96IY4
GeneID	1361;
Calculated MW	47kDa
Target Sequence	FQSGQVLAALPRTSRQVQLQNLTTTYEIVLWQPVTADLIVKKKQVHFFVNASDVDNVKAHLNVSGIPCSVLLA DVEDLIQQISNDTVSPRASASYEQYHSLNEIYSWIEFITERHPDMLTKIHIGSSFKEYPLYVLKVSQKEQAAK NAIWIDCGIHAREWISPAFCLWFIGHTQFYGIIGQYTNLLRLVDFYVMPVVNVGDYDYSWKKNRMWRKNRSFY ANNHCIGTDLNRNFASKHWCEEGASSSSCSETYCYGLYPESEPEVKAVASFLRRNINQIKAYISMHSYSQHIVFP YSYTRSKSKDHEELSLVASEAVRAIEKTSKNTRYTHGHGSETLYLAPGGGDDWIYDLGIKYSFTIELRDTGTYG FLLPERYIKPTCREAFAAVSKIAWHVIRNVLDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,1mM ZnCl ₂ ,10%Glycerol,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only