

# Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase FKBP7,FKBP7 (C-6His)



Catalog No: #AP75458

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75458-1 10ug #AP75458-2 50ug #AP75458-3 500ug #AP75458-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Peptidyl-Prolyl Cis-Trans Isomerase FKBP7,FKBP7 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9Y680
Uniprot	Q9Y680
GeneID	51661;
Calculated MW	23.94kDa
Target Sequence	QRQKKEESTEEVKIEVLRPENCSKTSKKGDLLNAHYDGYLAKDGSKFYCSRTQNEGHPKWFVLGVGQVIKG LDIAMTDMCPGKRRKVVIPPSFAYGKEGYAEGKIPDATLIFEIELYAVTKGPRSIETFKQIDMDNDRQLSKAEIN LYLQREFEKDEKPRDKSYQDAVLEDIFKKNHDHGDGFISPKEYNVYQHDELVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM GaCl <sub>2</sub> , 10%Glycerol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only