

Recombinant Human Tripeptidyl-Peptidase I, TPP1 (C-6His)

Catalog No: #AP75461

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75461-1 10ug #AP75461-2 50ug #AP75461-3 500ug #AP75461-4 1mg

Description

Product Name	Recombinant Human Tripeptidyl-Peptidase I, TPP1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	O14773
Uniprot	O14773
GeneID	1200;
Calculated MW	60.35kDa
Target Sequence	<p>SYSPEPDQRRTLPPGWVSLGRADPEEELSLTFALRQQNVERLSELVQAVSDPSSPQYGKYLLENVADLVRP SPLTLHTVQKWLLAAGAQAQKCHSVITQDFLTCWLSIRQAELLPLGAEFHVYVGGPTETHVVRSPHPYQLPQALA PHVDFVGGHLHRFPPTSSLRQRPEPQVTGTVGLHLGVTPSVIRKRYNLTSQDVGSGTSNNSQACAQFLEQYFH DSDLAQFMRLFGGNFAHQASVARVVGQQGRGRAGIEASLDVQYLMSAGANISTWVYSSPGRHEGQEPFLQ WLMLLSNESALPHVHTVSYGDDEDSLSSAYIQRVNTELMKAARGLTLLFASGDSGAGCWSVSGRHEFRPTF PASSPYVTTVGGTSFQEPFLITNEIVDYISGGGFSNVFPRPSYQEEAVTKFLSSSPHLPSSYFNASGRAYPDV AALSDGYWVVSNRVPIPWVSGTSASTPVFGGILSLINEHRILSGRPPLGFLNPRLYQQHGAGLFDVTRGCHES CLDEEVEGQGFCSGPGWDPVTGWGTPNFPALLKTLNLPVDHHHHHH</p>
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM GaCl ₂ , 10%Glycerol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only