

Recombinant Human Lysosomal Pro-X Carboxypeptidase,PRCP (C-6His)

Catalog No: #AP75462

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75462-1 10ug #AP75462-2 50ug #AP75462-3 500ug #AP75462-4 1mg

Description

Product Name	Recombinant Human Lysosomal Pro-X Carboxypeptidase,PRCP (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P42785
Uniprot	P42785
GeneID	5547;
Calculated MW	54.5kDa
Target Sequence	LRPALRALGSLHLPTNPTSLPAVAKNYSVLYFQQKVDHFGFNTVKTFFNQRYLVADKYWKKNGGSILFYTGNEG DIIWFCNNTGFMWDVAEELKAMLVFAEHRYYGESLPFGDNSFKDSRHLNFLTSEQALADFAELIKHLKRTIPGA ENQPVIAIGGSYGGMLAAWFRMKYPHMVVGALAASAPIWQFEDLVPCGVFMKIVTTDFRKSGPHCSSEIHR WDAINRLSNTGSLQWLTGALHLCSP LTSQDIQHLKDWISETWVNLAMVDYPYASNFLQPLPAWPIKVVCQYL KNPNVSDSLLLQNIQALNVYNYSGQVKCLNISETATSSLGTLGWSYQACTEVVMPFCTNGVDDMFEPHSW NLKELSDDCFQQWGVRRPSPWITMYGGKNISSHTNIVFNGELDPWSSGGVTKDITDTLVAVTISEGAHHL LRTKNALDPMSVLLARSLEVRHMKNWIRDFYDSAGKQHVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM NaAc-HAc, 150mM NaCl, 10% Glycerol, pH4.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only