

Recombinant Human Elafin, Trappin-2 (C-6His)

Catalog No: #AP75464



Package Size: #AP75464-1 10ug #AP75464-2 50ug #AP75464-3 500ug #AP75464-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Elafin, Trappin-2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P19957
Uniprot	P19957
GeneID	5266;
Calculated MW	11kDa
Target Sequence	AVTGVPVKGQDTVKG RVPFNGQDPVKGQVSVKGQDKVKAQEPVKG PVSTKPGSCPIILIRCAMLNPPNRCLK DTDCPGIKKCCGEGSCGMACFVPQVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, pH 7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only