

Recombinant Human ERO1-Like Protein α , ERO1L (C-6His)

Catalog No: #AP75466

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75466-1 10ug #AP75466-2 50ug #AP75466-3 500ug #AP75466-4 1mg

Description

Product Name	Recombinant Human ERO1-Like Protein α , ERO1L (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96HE7
Uniprot	Q96HE7
GeneID	30001;
Calculated MW	53kDa
Target Sequence	EEQPPETAAQRCFCQVSGYLDDCTCDVETIDRFNNYRLFPRQLKLLSESDYFRYYKVNLRPCPFWNDISQCG RRDCAVKPCQSDEVPDGIKSASYKYSEEANLIEECEQAERLGAVDESLSEETQKAVLQWTKHDDSSDNFCE ADDIQSPEAEYVDLLNPERYTYGKGPDAWKIWNVIYEENCFKPQTIKRPLNPLASGQGTSEENTFYSWLEGL CVEKRAFYLISGLHASINVHLSARYLLQETWLEKKWGHNITEFQQRFDGILTEGEGPRRLKNLYFLYLIELRAL SKVLPFFERPDFQLFTGNKIQDEENKMLLLEILHEIKSFPLHFDENSFFAGDKKEAHKLKEDFRLHFRNISRIMD CVGCFKCRLWGKLQTQGLGTALKILFSEKLIANMPESGPSYEFHLTRQEIVSLFNAFGRISTSVKELENFRNLL QNIHVDHHHHHH
Formulation	Lyophilized from a 0.2 μ m filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only