

Recombinant Human Reticulocalbin-3,RCN3 (C-6His)

Catalog No: #AP75468



Package Size: #AP75468-1 10ug #AP75468-2 50ug #AP75468-3 500ug #AP75468-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Reticulocalbin-3,RCN3 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96D15
Uniprot	Q96D15
GeneID	57333;
Calculated MW	36.2kDa
Target Sequence	KPSPDAGPHGQGRVHQAAPLSAPHDDAHGNFQYDHEAFLGREVAKEFDQLTPEESQARLGRIVDRMDRA GDGDGWVSLAELRAWIAHTQQRHIRDSVSAAWDTYDTRDGRVGEELRNATYGHYAPGEEFHDVEDAET YKMLARDERRFRVADQDGDSDMATREELTAFLHPEEFPHMRDIVIAETLEDLDRNKDGYVQVEEYIADLYSAE PGEEPAWVQTERQQFRDLNKGHLDGSEVGHWWLPPAQDQPLVEANHLLHESDTRDKDGRLSKAEILG NWNMFVGSQATNYGEDLTRHHDELVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 1mM DTT, 10%Glycerol, pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only