

# Recombinant Human Cadherin-11,OB-Cadherin,CDH11 (C-6His)

Catalog No: #AP75469

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75469-1 10ug #AP75469-2 50ug #AP75469-3 500ug #AP75469-4 1mg

## Description

Product Name	Recombinant Human Cadherin-11,OB-Cadherin,CDH11 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96CZ9
Uniprot	Q96CZ9
Calculated MW	66.63kDa
Target Sequence	FAPERRGHLRPSFHGHHEKGGKQVLQRSKRGWVWNQFFVIEEYTGPDVPLVGRGLHSDIDSGDGNIKYILSG EGAGTIFVIDDKSGNIHATKTLDRREERAQYTLMAQAVDRDTRNRPLEPPSEFIVKVQDINDNPPEFLHETYHANV PERSNVGTSVIQVTASDADDPTYGNSAKLVYSILEGQPYFSVEAQTGIIRTALPNMDREAKEEYHVVIQAKDMG GHMGGLSGTTKVMITLTDVNDNPPKFPQSVYQMSVSEAAVPGEEVGRVKAKDPDIGENGLVTYNIVDGDGME SFEITTDYETQEGVIKLLKPVDFETKRAYSLKVEAANVHIDPKFISNGPFKDTVTVKIAVEDADEPPMFLAPSYIH EVQENAAAGTVVGRVHAKDPDAANSPIRYSIDRHTDLDRFFINPEDGFIKTKPLDREETAWLNITVFAAEIHN RHQEAKVPVAIRVLDVNDNAPKFAAPYEGFICESDQTKPLSNQPIVTISADDKDDTANGPRFIFSLPPEIIHNPNF TVRDNRDNTAGVYARRGGFSRQKQDLYLLPIVISDGGIPMSSTNTLTIKVCGCDVNGALLSCNAEAYILNAGL STVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only