

Recombinant Human Leukocyte Ig-Like Receptor A2,LILRA2,ILT1,CD85h (C-6His)

Catalog No: #AP75470

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75470-1 10ug #AP75470-2 50ug #AP75470-3 500ug #AP75470-4 1mg

Description

Product Name	Recombinant Human Leukocyte Ig-Like Receptor A2,LILRA2,ILT1,CD85h (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8N149
Uniprot	Q8N149
GeneID	11027;
Calculated MW	44.8kDa
Target Sequence	GHLPKPTLWAEPGSVIIQGSPVTLRCQGSLQAEYHLYRENKSASWVRRIQEPGKNGQFPPIPSITWEHAGRYH CQYYSHNHSSEYSDPLELVVTGAYSKPTLSALPSPVVTLGGNVTLCQVSVQVAFDGFILCKEGEDHPQRLNSH SHARGWSWAIFSVGPVSPSRRWSYRCYAYDSNSPYVWSLPSDLELLVPGVSKKPSLSVQPGPMVAPGESL TLQCVSDVGYDRFVLYKEGERDFLQRPGWQPQAGLSQANFTLGPVSPSHGGQYRCYSAHNLSSEWSAPSD PLDILITGQFYDRPSSLVQPPTVAPGKNVTLLCQSRGQFHTFLLTKEGAGHPPLHLRSEHQAAQNQAEFRMG PV TSAHVGTYRCYSSLSNPYLLSLPSDPLELVVSASVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only