

# Recombinant Human Prostaglandin-H2 D-Isomerase,PTGDS (C-6His)

Catalog No: #AP75473

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75473-1 10ug #AP75473-2 50ug #AP75473-3 500ug #AP75473-4 1mg

## Description

Product Name	Recombinant Human Prostaglandin-H2 D-Isomerase,PTGDS (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P41222
Uniprot	P41222
GeneID	5730;
Calculated MW	19.7kDa
Target Sequence	APEAQVSVQPNFQQDKFLGRWFWSAGLASNSSWLREKKAALSMCKSVVAPATDGGLNLTSTFLRKNQCETRT MLLQPAGSLGSSYSRSPHWGSTYSVSVVETDYDQYALLYSQGSKGPGEDFRMATLYSRTQTPRAELKEKFT AFCKAQGFTEDTIVFLPQTDKCMTEQVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl,150mM NaCl,10%Glycerol,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only