

# Recombinant Human Brevican Core Protein, BCAN, BEHAB (C-6His)



Catalog No: #AP75487

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75487-1 10ug #AP75487-2 50ug #AP75487-3 500ug #AP75487-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Brevican Core Protein, BCAN, BEHAB (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q96GW7
Uniprot	Q96GW7
GeneID	63827;
Calculated MW	97.75kDa
Target Sequence	DVLEGDSSDRAFRVRIAGDAPLQGVLGALTIPCHVHYLRPPPSRRAVLGSPRVKWTFLSRGREAEVLVAR GVRVKVNEAYRFRVALPAYPASLTDVSLALSELRPNDSGIYRCEVQHIGDSSDAVEVKVKGWVFLYREGSAR YAFSFGAQEACARIGAHIAITPEQLYAAYLGGYEQC DAGWLS DQTVRYPIQTPREACYGDMDFGVRNYG VVDPDDLVDVYCYAEDLNGELFLGDPPEKLTLEEARAYCQERGAEIATTGQLYAAWDGGLDHCSPGWLADG SVRYPIVTPSQRCCGGGLPGVKTLFLFPNQTGFPNKHSRNFVYCFRDSAQPSAIP EASN PASNDG LEAIVT VTETLEELQLPQEATESESRGAIYSIPIMEDGGGGSSSTPEDPAEAPRTLLEFETQSMVPPTGFSEEEGKALEEE EKYEDEEEKEEEEEVEDEALWAWPSELSSPGPEASLPTEPAAQEKSLSQAPARAVLQPGASPLPDGESE ASRPPRVHGPPTETLPTPRERNLASPSTLVEAREVGEATGGPELSGVPRGESEETGSSEGAPSLPATRA PEGTRELEAPSEDNSGRTAPAGTSVQAQPVLPTDSASRGGVAVVPASGDCVPSCHNGGTCL EEEEGVRCL CLPGYGGDLCDVGLRFCNPGWDAFQGACYKHFSTRRSWEEAETQCRMYGAHLASISTP EEQDFINNRYREY QWIGLNDRTIEGDFLWSDGVPLLYENWNPGQPDSYFLSGENCVVMVWHDQQWSDVPCNYHLSYTCCKMGL VSCGPPPELPLAQVFGRPRLRYEVDTVLRYRCREGLAQRNLPLIRCQENGRWEAPQISCVPRRPARALHPEE DPEGRQGRLLGRWKALLIPPSSPMPGPHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only