

Recombinant Human Activin Receptor 1A,Activin RI,ALK-2,ACVR1 (C-6His)

Catalog No: #AP75489

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75489-1 10ug #AP75489-2 50ug #AP75489-3 500ug #AP75489-4 1mg

Description

Product Name	Recombinant Human Activin Receptor 1A,Activin RI,ALK-2,ACVR1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q04771
Uniprot	Q04771
GeneID	90;
Calculated MW	12.6kDa
Target Sequence	MEDEKPKVNPPLYMCVCEGLSCGNEDHCEGQQCFSSLSINDGFHVYQKGCQVVEEQGKMTCKTPSPGQA VECCQGDWCNRNITAQLPTKGKSFPGTQNFHLEVVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only