

Recombinant Human Neuroligin 4, X-Linked,NLGN4X (C-6His)



Catalog No: #AP75490

Orders: order@signalwayantibody.com

Package Size: #AP75490-1 10ug #AP75490-2 50ug #AP75490-3 500ug #AP75490-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Neuroligin 4, X-Linked,NLGN4X (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8N0W4
Uniprot	Q8N0W4
GeneID	57502;
Calculated MW	72.46kDa
Target Sequence	QAQYPVVNTNYGKIRGLRTPNEILGPVEQYLGVPYASPTGERRFQPPEPPSSWTGIRNTTQFAAVCPQHL DERSLLHDMLPIWFTANLDTLMTYVQDQNEDECLYLNIVPTEDDIHQNSKKPVMVYIHGGSYMEGTGNMIDG SILASYGNVIVITINYRLGILGFLSTGDQAAKGYGLLDQIQALRWIEENVGAFGGDPKRVTFGSGAGASCVSLL TLSHYSEGLFQKAIQSGTALSSWAVNYQPAKYTRILADKVGCMMLDTTDMVECLRKNKYKELIQQTITPATYHI AFGPVIDGDVIPDDPQILMEQGEFLNYDIMLGVNQGEGLKFVDGIVDNEGDVTPNDFDFSVSNFVDNLYGYPE GKDTLRETIKFMYPDWADKENPETRRKTLVALFTDHQWVAPAVATADLHAQYGSPTYFYAFYHHCQSEMKPS WADSAHGDEVYVFGIPMIGPTELFSCNFSKNDVMLSAAVVMTYWTNFAKTGDPNQPVQDTKFIHTKPNRFE EVAWSKYNPKDQLYLHIGLKPRVRDHYRATKVAFWLELPHLHNLNEIFQYVSTTTKVPPDMSFPYGTRRS PAKIWPTRKPAITPANNPKHSDPHKTPEDTTVLIETKRDYSTEIVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM PB, 150mM NaCl, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks. Reconstituted protein solution can be stored at 4-7C for 2-7 days. Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only