

Recombinant Human Vitamin K-Dependent Protein C,PROC (C-6His)



Catalog No: #AP75491

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75491-1 10ug #AP75491-2 50ug #AP75491-3 500ug #AP75491-4 1mg

Description

Product Name	Recombinant Human Vitamin K-Dependent Protein C,PROC (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P04070
Uniprot	P04070
GeneID	5624;
Calculated MW	51.1kDa
Target Sequence	TPAPLDSVFSSSERAHQVLRIKRRANSFLEELRHSSLERECIEEICDFEEAKEIFQNVDDTLAFWSKHVDGDQC LVLPLEHPCASLCCGHGTCIDGIGSFSCDCRSWEGRFQREVSVFLNCSLDNGGCTHYCLEEVGWRRCSA PGYKLGDDLLQCHPAVKFPCGRPWKRMKKRSHLKRDTEDQEDQVDPRLIDGKMTRRGDSPWQVLLDSK KKLACGAVLIHPSWVLTAAHCMDESKLLVRLGEYDLRRWEKWELDLDIKEVVFHPNYSKSTTDNDIALHLA QPATLSQTIVPICLPDSGLAERELNQAQGETLVGTWGYHSSREKEAKRNRTFVLFNIKIPVPHNECSEVMSN MVSENMLCAGILGDRQDACEGDSGGPMVASFHGTWFLVGLVSWGEGCGLLHNYGVYTKVSRYLWDWIHGHIR DKEAPQKSWAPLEHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10%Glycerol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only