

# Recombinant Human Hyaluronidase-1, HYAL1 (C-6His)

Catalog No: #AP75492

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75492-1 10ug #AP75492-2 50ug #AP75492-3 500ug #AP75492-4 1mg

## Description

Product Name	Recombinant Human Hyaluronidase-1, HYAL1 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q12794
Uniprot	Q12794
GeneID	3373;
Calculated MW	47.18kDa
Target Sequence	FRGPLLPNRPFTTVWNANTQWCLERHGVDDVSVFDVVANPGQTFRGPDMTIFYSSQLGTYPYYTPTGEPVF GGLPQNASLIAHLARTFQDILAAIPAPDFSGLAVIDWEAWRPRWAFNWDTKDIYRQRSRALVQAQHPDWPAP QVEAVAQDQFQGAARAWMAGTLQLGRALRPRGLWGFYGFPCYNYDFLSPNYTGQCPSGIRAQNDQLGWL WGQSRALYPSIYMPAVLEGTKSQMYVQHRVAEAFRVAVAAGDPNLPVLPYVQIFYDTTNHFLPLDELEHSL GESAAQGAAGVVLWVSWENTRTKESCQAIKEYMDTTLGPFILNVTSGALLCSQALCSGHGRCVVRTSHPKAL LLLNPAFSIQLTPGGGPLSLRGALSLEDQAQMAVEFKCRCYPGWQAPWCERKSMWVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM TrisHCl, 150mM NaCl, 10%Glycerol, pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only