

# Recombinant Human Fms-Related Tyrosine Kinase 3 Ligand,FLT3LG (C-6His)

Catalog No: #AP75493

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75493-1 10ug #AP75493-2 50ug #AP75493-3 500ug #AP75493-4 1mg

## Description

Product Name	Recombinant Human Fms-Related Tyrosine Kinase 3 Ligand,FLT3LG (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P49771
Uniprot	P49771
GeneID	2323;
Calculated MW	18.9kDa
Target Sequence	TQDCSFQHSPISDFAVKIRELSDYLLQDYPVTVASNLQDEELCGGLWRLVLAQRWMERLKTVAGSKMQGLL ERVNTEIHFVTKCAFQPPPSCLRFVQTNISRLLETSEQLVALKPWITRQNFSRCLELQCQPDSSTLPPPWSP RPLEATAPTAPQPVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of 20mM TrisHClpH8.0.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only