

# Recombinant Human Proprotein Convertase Subtilisin,Kexin Type 9,PCSK9 (C-AVI)



Catalog No: #AP75495

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75495-1 10ug #AP75495-2 50ug #AP75495-3 500ug #AP75495-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Proprotein Convertase Subtilisin,Kexin Type 9,PCSK9 (C-AVI)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q8NBP7
Uniprot	Q8NBP7
GeneID	255738;
Calculated MW	71.05kDa
Target Sequence	QEDEDGDYEELVLALRSEEDGLAEAPEHGTTATFHRCADPWRLPGTYVVVLKEETHLSQSERTARRLQAQA ARRGYLTKILHVFHGLLPGLVKMSGDLLELALKLPHVDYIEEDSSVFAQSIPWNLERITPPRYRADEYQPPDG GSLVEVYLLDTSIQSDHREIEGRVMVTDFFENVPEEDGTRFHRQASKCDSHGTHLAGVVSGRDAGVAKGASM RSLRVLNCQGGKGTVSGTLIGLEFIRKSQLVQPVGPLVLLPLAGGYSRVLNAACQRLARAGVWLTAAGNFRD DACLYSPASAPEVITVGATNAQDQPVTGLTGTNFGRCVDLFPAGEDIIIGASSDCSTCFVSQSGTSQAAAHVA GIAAMMLSAEPELTLAELRQLIHFSKDVINEAWFPEDQRVLTPNLVAALPPSTHGAGWQLFCRTVWSAHSG PTRMATAIARCAPDEELLSCSSFSRSGKRRGERMEAQGGKLVCRAHNAFGGEGVYAIARCCLLPQANCSVHT APPAEASMGTRVHCHQGGHVLGCSHWEVEDLGTHKPPVLRPRGQPNQCVGHREASIHASCCHAPGLEC KVKEHGIPAPQEQVTVACEEGWTLTGCSALPGTSHVLGAYAVDNTCVVRSRDVSTTGSTSEEAVTAVAICCR SRHLAQASQELQEFGGGLNDIFEAQKIEWHE
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only