

Recombinant Mouse Proprotein Convertase Subtilisin,Kexin Type 9,PCSK9 (C-6His)



Catalog No: #AP75496

Orders: order@signalwayantibody.com

Package Size: #AP75496-1 10ug #AP75496-2 50ug #AP75496-3 500ug #AP75496-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Mouse Proprotein Convertase Subtilisin,Kexin Type 9,PCSK9 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q80W65
Uniprot	Q80W65
GeneID	100102;
Calculated MW	74.8kDa
Target Sequence	<p>QDEDGDYEELMLALPSQEDGLADEAAHVATATFRRCsKEAWRLPGTYIVVLMETQRLQIEQTAHRLQTRAA RRGYVIKVLHIFYDLFPGFVKMSSDLLGLALKLPHVEYIEEDSFVFAQSIPWNLERIIPAWHQTEEDRSPDGSS QVEVYLLDTSIQGAHREIEGRVTITDFNSVPEEDGTRFHRQASKCDShGTHLAGVVSGRDAGVAKGTSLSLRL VLNCQKGTVSGTLIGLEFIRKSQLIQPSGPLVLLPLAGGYSRILNAACRHLARTGVVLAAGNFRDDACLY SPASAPEVITVGATNAQDQPVTLGTLGTNFGRCVDFAPGKDIIGASSDCSTCFMSQSGTSQAAAHVAGIVAR MLSREPTLTLAELRQRLIHfSTKdVINMAWFPEDQQVLTpNLVATLPPSTHETGGQLLCRTVWSAHSGPTRTA TATARCAPEEELLSCSSFSRSGRRRGDWIEAIGGQQVCKALNAFGGEGVYAVARCCLVPRANCSIHNTPAAR AGLETHVHCHQKDHVLTGCSFHWEVEDLSVRRQPALRSRRQPGQCVGHQAASVYASCCHAPGLECKIKEH GISGPSEQVTVACEAGWTLTGCNVLPGASLTlGAYSVDNLCVARVHDTARADRTSGEATVAAAICCRSRPSA KASWVQHSHHHHH</p>
Formulation	Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, pH 7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only