

Recombinant *Macaca nemestrina* PCSK9(C-6His)

Catalog No: #AP75497



Package Size: #AP75497-1 10ug #AP75497-2 50ug #AP75497-3 500ug #AP75497-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

| | |
|-----------------|--|
| Product Name | Recombinant <i>Macaca nemestrina</i> PCSK9(C-6His) |
| Host Species | Human Cells |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | <i>Macaca nemestrina</i> |
| Accession No. | A8T662 |
| Uniprot | A8T662 |
| Calculated MW | 14&59kDa |
| Target Sequence | <p>QEDEDGDYEELVLALRSEEDGLADAPEHGATATFHRCADPWRLPGTYVVVLKEETHRSQSERTARRLQAQ AARRGYLTKILHVFHLLPGFLVKMSGDLLELALKLPHVDYIEEDSSVFAQSIPWNLERITPARYRADEYQPPKG GSLVEVYLLDTSIQSDHREIEGRVMVTDVFESVPEEDGTRFHRQASKCDSHGTHLAGVVSGRDAGVAKGAGLR SLRVLNCQGGKGTVSGTLIGLEFIRKSQLVQPVGPLVLLPLAGGYSRVFNAACQLRARAGVVLVTAAGNFRDD ACLYSPASAPEVITVGATNAQDQPVTGLTGLTNFGRCVDLFAPGEDIIGASSDCSTCFVSRSGTSQAAAHVAGI AAMMLSAEPELTLAELRQLIHFSKDVINEAWFPEDQRVLTPNLVAALPPSTHRAGWQLFCRTVWSAHS GP TRMATAVARCAQDEELLSCSSF SRSGKRRGERIEAQGGKRV CRAHNAFGGEGVYAIARCCLLPQVNC SVHT APPAGASMGTRVHCHQQGHVLTGCSHWEEVDLGT HKPPVLRPRGQPNQCVGHREASIHASCCHAPGLEC KVKEHGIPAPQEQVIVACEDGWTLTGCSALPGTSHVLGAYAVDNTCVVRSRDVSTTGSTSEEVAVAVAICCRS RHLVQASQELQH HHHHHH</p> |
| Formulation | Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 20% Glycerol, pH 7.4. |
| Storage | Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only