

Recombinant Human Proprotein Convertase Subtilisin,Kexin Type 9,PCSK9 (D374Y, C-6His)



Catalog No: #AP75499

Orders: order@signalwayantibody.com

Package Size: #AP75499-1 10ug #AP75499-2 50ug #AP75499-3 500ug #AP75499-4 1mg

Support: tech@signalwayantibody.com

Description

| | |
|-----------------|--|
| Product Name | Recombinant Human Proprotein Convertase Subtilisin,Kexin Type 9,PCSK9 (D374Y, C-6His) |
| Host Species | Human Cells |
| Purification | Greater than 95% as determined by reducing SDS-PAGE. |
| Specificity | Human |
| Accession No. | Q8NBP7 |
| Uniprot | Q8NBP7 |
| GeneID | 255738; |
| Calculated MW | 71.1kDa |
| Target Sequence | QEDEDGDYEELVLALRSEEDGLAEPEHGTTATFHRCADPWRLPGTYVVVLKEETHLSQSERTARRLQAQA ARRGYLTKILHVFHGLLPGLVKMSGDLLELALKLPHVDYIEEDSSVFAQSIPWNLERITPPRYRADEYQPPDG GSLVEVYLLDTSIQSDHREIEGRVMVTDFFENVPEEDGTRFHRQASKCDSHGTHLAGVVSGRDAGVAKGASM RSLRVLNCQGGKGTVSGTLIGLEFIRKSQLVQPVGPLVLLPLAGGYSRVLNAACQRLARAGVWLTAAGNFRD DACLYSPASAPEVITVGATNAQDQPVTGLTGTNFGRCVDLFPAGEDIIIGASSYCSTCFVSQSGTSQAAAHVA GIAAMMLSAEPELTLAELRQLIHFSKDVINEAWFPEDQRVLTPNLVAALPPSTHGAGWQLFCRTVWSAHSG PTRMATAIARCAPDEELLSCSSFSRSGKRRGERMEAQGGKLVCRAHNAFGGEGVYAIARCCLLPQANCSVHT APPAEASMGTRVHCHQGGHVLGCSHWEVEDLGTHKPPVLRPRGQPNQCVGHREASIHASCCHAPGLEC KVKEHGIPAPQEQVTVACEEGWTLTGCSALPGTSHVLGAYAVDNTCVVRSRDVSTTGSTSEEAVTAVAICCR SRHLAQASQELQHSHHHHH |
| Formulation | Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, pH 7.4. |
| Storage | Store at < -20C, stable for 6 months after receipt. Please minimize freeze-thaw cycles. |

Note: This product is for in vitro research use only