

Recombinant Human IL-1 Receptor Accessory Protein,IL-1RAcP,IL-1R3 (C-Fc-6His)



Catalog No: #AP75500

Orders: order@signalwayantibody.com

Package Size: #AP75500-1 10ug #AP75500-2 50ug #AP75500-3 500ug #AP75500-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human IL-1 Receptor Accessory Protein,IL-1RAcP,IL-1R3 (C-Fc-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9NPH3
Uniprot	Q9NPH3
GeneID	3556;
Calculated MW	39.7kDa
Target Sequence	SERCDDWGLDTMRQIQVFEDPARIKCPLFEHFLKFNYSTAHSAGLTLIWIYWTRQDRDLEEPINFRLPENRISK EKDVLWFRPTLLNDTGNVTCMLRNNTTYCSKVAFFLEVVQKDFNSPMKLPVHKLYIEYGIQRITCPNVDGYF PSSVKPTITWYMGCYKIQNFNNVIPEGMNLFLIALISNNGNYTCVVVYPENGRTFHLTRTLTVKVVGSPKNAV PPVIHSPNDHVVEKEPGEELLIPCTVYFSFLMDSRNEVWWTIDGKKPDDITIDVTINESISHSRTEDETRTQILS IKKVTSEDLKRSYVCHARSAKGEVAKAAKVKQKGNRCGQVDDIEGRMDEPKSCDKTHTCPPCPAPPELLGGPS VFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQ DWLNGKEYCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESN GQPENNYKTTTPVLDSGFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQKSLSLSPGKHHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only