

Recombinant Human IL-20 Receptor Subunit α

α ,IL20RA (C-Fc)

Catalog No: #AP75501

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75501-1 10ug #AP75501-2 50ug #AP75501-3 500ug #AP75501-4 1mg

Description

Product Name	Recombinant Human IL-20 Receptor Subunit α ,IL20RA (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q9UHF4
Uniprot	Q9UHF4
GeneID	53832;
Calculated MW	26.3kDa
Target Sequence	VPCVSGGLPKPANITFLSINMKNVLQWTPPEGLQGKVKVITYTVQYFIYGQKKWLNKSECRNINRTYCDLSAETS DYEHQYYAKVKAIWGTKCSKWAESGRFYPFLETQIGPPEVALTTDEKSISVVLTAPEKWKRNPEDLPVSMQQI YSNLKYNVSVLNTKSNRTWSQCVTNHTLVLTWLEPNTLYCVHVESFVPGPPRAQPSEKQCARTLKDQSSEF KAKVDDIEGRMDEPKSCDKTHTCPPCPAPPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVDVSHEDPEVKFN WYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQV YTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSGDSFFLYSKLTVDKSRWQQG NVFSCSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 μ m filtered solution of PBS, pH 7.4 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only