

Recombinant Human NKG2-A,NKG2-B Type II Integral Membrane Protein(N-8His)

Catalog No: #AP75502

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75502-1 10ug #AP75502-2 50ug #AP75502-3 500ug #AP75502-4 1mg

Description

Product Name	Recombinant Human NKG2-A,NKG2-B Type II Integral Membrane Protein(N-8His)
Host Species	Human cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P26715
Uniprot	P26715
GeneID	3821;
Calculated MW	16.5kDa
Target Sequence	HHHHHHHRHNSSLNTRTQKARHCGHCPEEWITYSNSCYYIGKERRTWEESSLACTSKNSSLLSIDNEEEM KFLSIISPSSWIGVFRNSSHHPWVTMNGLAFKHEIKDSDNAELNCAVLQVNRLKSAQCGSSIIYHCKHKL
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only