

# Recombinant Human Pregnancy-Specific β-1-Glycoprotein 9,PSBG9 (C-6His)



Catalog No: #AP75503

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75503-1 10ug #AP75503-2 50ug #AP75503-3 500ug #AP75503-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human Pregnancy-Specific β-1-Glycoprotein 9,PSBG9 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q00887
Uniprot	Q00887
GeneID	5678;
Calculated MW	45.6kDa
Target Sequence	EVTIEAQPPKVSEGKDVLLLVHNLPGYFWYKGMTDLHYHYIISYIVDGKIIIGPAYSGRETIVYSNASLLIQ NVTRKDAGTYTLHIIKRGDETREEIRHFTFTLYLETPKPYISSNLNPREAMEAVRLICDPETLDASYLWWMNG QSLPVTHRLQLSKTNRTLYLFGVTKYIAGPYECEIRNPVSASRSDPVTLLNLLPKLPIPYITINLNPRENKDVLAF TCEPKSENYTYIWWLNGQSLPVSPGVKRIENRILIPSVTRNETGPYQCEIQDRYGGLRSNPVILNVLYGPD PRIYPSFTYYRSGENLDLSCFTESNPPAEYFWTINGKFQQSGQKLFIPQITRNHSGLYACSVHNSATGKEISKS MTVKVSGPCHGDLTESQSVDHHHHHH
Formulation	Supplied as a 0.2 μm filtered solution of 20mM PB,150mM NaCl,pH7.5.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only