

# Recombinant Human N-Acetylgalactosamine kinase,GALK2 (C-6His)

Catalog No: #AP75505

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Package Size: #AP75505-1 10ug #AP75505-2 50ug #AP75505-3 500ug #AP75505-4 1mg

## Description

Product Name	Recombinant Human N-Acetylgalactosamine kinase,GALK2 (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q01415
Uniprot	Q01415
GeneID	2585;
Calculated MW	51.4kDa
Target Sequence	MATESPATRRVQVAEHPRLKEMFNKFGSIPKFYVRAPGRVNIIGEHIDYCGYSVLPMAVEQDVLIAVEPV KTYALQLANTNPLYPDFSTSANNIQIDKTKPLWHNYFLCGLKGIQEHFGLSNLTGMNCLVDGNIPPSSGLSSSS ALVCCAGLVTLTVLGRNLSKVELAEICAKSERYIGTEGGGMDQSSISFLAEEGTAKLIEFSPLRATDVKLPSGAVF VIANSCVEMNKAATSHFNIRVMECRLAAKLLAKYKSLQWQDKVLRLEEVQAKLGLISLEEMLLVTEALHPEPYNP EEICRCLGISLEELRTQILSPNTQDVLIFKLYQRAKHVYSEAAARVLQFKKICEEAPENMVQLLGELMNQSHMSC RDMYECSCPELDQLVDICRKFQAQGSRLTGAGWGGCTVSMVPADKLPFLANVHKAYYQRS DGLAPEKQS LFATKPGGGALVLEAVDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of 20mM PB,150mM NaCl,pH7.4.
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only