

Recombinant Human Alkaline Phosphatase,ALPP,PLAP (C-6His)



Catalog No: #AP75508

Orders: order@signalwayantibody.com

Package Size: #AP75508-1 10ug #AP75508-2 50ug #AP75508-3 500ug #AP75508-4 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Alkaline Phosphatase,ALPP,PLAP (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P05187
Uniprot	P05187
GeneID	250;
Calculated MW	53.7kDa
Target Sequence	IIPVEEENPDFWNREAAEALGAAKKLQPAQTAANKLIIFLGDGMGVSTVTAARILKGQKKDKLGP EIPLAMDRFP YVALSKTYNVDKHVPDSGATATAYLCGVKGNFQTIGLSAAARFNQCNTTRGNEVISVMNRAKKAGKSVGVVT TTRVQHASPAGTYAHTVNRNWYSDADV PASARQEGCQDIATQLISNMDIDVILGGGRKYMFRMGTPDPEYPD DYSQGGTRLDGKNLVQEWLAKRQGARYVWNRT ELMQASLDPSVTHLMGLFEPGDMKYEIHRDSTLDPSLME MTEAALRLLSRNPRGFFLFVEGGRIDHGHESRAYRALTETIMFDDAIERAGQLTSEEDT LSLVTADHSHVFSF GGYPLRGSSIFGLAPGKARDRKAYTVLLYGNGPGYVLKDGARPDVTESESGSPEYRQQSAVPLDEETHAGED VAVFARGPQAHLVHGVQEQT FIAHVMAFAACLEPYTACDLAPPAGTTD VDHHHHHH
Formulation	Supplied as a 0.2 µm filtered solution of PBS, pH 7.4 .
Storage	Store at < -20C, stable for 6 months after receipt.Please minimize freeze-thaw cycles.

Note: This product is for in vitro research use only