

Recombinant Human Lymphocyte Function-associated Antigen 3,CD58,LFA-3 (C-Fc)

Catalog No: #AP75509

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75509-1 10ug #AP75509-2 50ug #AP75509-3 500ug #AP75509-4 1mg

Description

Product Name	Recombinant Human Lymphocyte Function-associated Antigen 3,CD58,LFA-3 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P19256
Uniprot	P19256
GeneID	965;
Calculated MW	48.6kDa
Target Sequence	FSQQIYGVVYGNVTFHVPSNVPLKEVLWKKQKDKVAELENSEFRAFSSFKNRVYLDTVSGSLTIYNLTSSDED EYEMESPNITDTMKFFLYVLES LPSPTLTCAL TNGSIEVQCM IPEYYNSHRGLIMYSWDCPMEQCKRNSTSIYF KMENDLPQKIQCTLSNPLFNTTSSIIL TTCIPSSGHSRHRVDDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSV FLFPPKPKDTLMISRTPEVTCVVVDVSHEDPEVKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQD WLNQKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNG QPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNV FSCSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only