

Recombinant Human CD7,Leu-9 (C-Fc)

Catalog No: #AP75510



Package Size: #AP75510-1 10ug #AP75510-2 50ug #AP75510-3 500ug #AP75510-4 1mg

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CD7,Leu-9 (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P09564
Uniprot	P09564
GeneID	924;
Calculated MW	43.4kDa
Target Sequence	AQEVQQSPHCTTVPVGASVNITCSTSGGLRGIYLRQLGPQPQDIYYEDGVVPTTDRRFRGRIDFSGSQDNLTI TMHRLQLSDTGTYTCAITEVNVYGSGLTLVLVTEEQSQGWHRCSAPPASALPAPPTGSALPDPQTASALP DPPAASALPDIEGRMDEPKSCDKTHTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVVVDVSHEDPE VKFNWYVDGVEVHNAKTKPREEQYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPR EPQVYTLPPSREEMTKNQVSLTCLVKGFYPSDIAVEWESNGQPENNYKTPPVLDSDGSFFLYSKLTVDKSR WQQGNVFSCSVMHEALHNHYTQKSLSLSPGK
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only