

Recombinant Human E-Selectin,,ELAM1,CD62E (C-6His)



Catalog No: #AP75511

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75511-1 10ug #AP75511-2 50ug #AP75511-3 500ug #AP75511-4 1mg

Description

Product Name	Recombinant Human E-Selectin,,ELAM1,CD62E (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	P16581
Uniprot	P16581
GeneID	6401;
Calculated MW	58.6kDa
Target Sequence	WSYNTSTEAMTYDEASAYCQQRythlVAIQNKEEIEYLNLSILSYSPSYWIGIRKVNNVVVWVGTQKPLTEEA KNWAPGEPNNRQKDEDCVEIYIKREKDVGMWNDERCSKKLALCYTAACTNTSCSGHGECVETINNYTCKC DPGFSGLKCEQIVNCTALESPEHGSLVCSHPLGNFSYNSSCSISCDRGLPSSMETMQCMSSGEWSAPIPAC NVVECDAVTNPANGFVECFQNPGSFPWNTTCTFDCEEGFELMGAQSLQCTSSGNWDNEKPTCKAVTCRAV RQPQNGSVRCSHPAGEFTFKSSCNFTCEEGLQGPQVVECTTQGGWTQQIPVCEAFQCTALSNPERGY MNCLPSASGSFRYGSSCEFSCEQGFVLKGSKRLQCGPTGEWDNEKPTCEAVRCDVAHQPPKGLVRCASPI GEFTYKSSCAFSCEEGLHGSTQLECTSQGQWTEEVPSQVVKCSSLAVPGKINMSCSGEPVFGTVCKFA CPEGWTLNGSAARTCGATGHWSGLLPTCEAPTESNIPVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, 4% Mannitol,pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only