

# Recombinant Mouse Granulocyte Colony-Stimulating Factor,G-CSF,CSF1(C-6His)



Catalog No: #AP75513

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #AP75513-1 10ug #AP75513-2 50ug #AP75513-3 500ug #AP75513-4 1mg

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Mouse Granulocyte Colony-Stimulating Factor,G-CSF,CSF1(C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	P09920
Uniprot	P09920
GeneID	12985;
Calculated MW	19.8kDa
Target Sequence	VPLVTVSALPPSLPLPRSFLLKSLEQVRKIQASGSVLLLEQLCATYKLCHPPEELVLLGHSLGIPKASLSGCSSQAL QQTQCLSQLHSGLCCLYQGLLQALSGISPALAPTLDLLQLDVANFATTIWQQMENLGVAPTQPTQSAMPAFTS AFQRRAGGVLAISYLQGFLETARLALHHLAHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH7.4.
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only