

Recombinant Human IL-10 Receptor Subunit β,IL-10RB (C-Fc)



Catalog No: #AP75547

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75547-1 10ug #AP75547-2 50ug #AP75547-3 500ug #AP75547-4 1mg

Description

Product Name	Recombinant Human IL-10 Receptor Subunit β,IL-10RB (C-Fc)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Human
Accession No.	Q08334
Uniprot	Q08334
GeneID	3588;
Calculated MW	50.6kDa
Target Sequence	MVPPPENVRMNSVNFKNILQWESPAFAKGNLFTTAQYLSYRIFQDKCMNTLTECDFSSLSKYGDHTLRVRAE FADEHSDWVNITFCPVDDTIIGPPGMQVEVLADSLHMRFLAPKIENEYETWTMKNVYNSWTYNVQYWKNGTD EKFKITPQYDFEVLRLNLEPWTTYCVQVRGFLPDRNKAGEWSEPVCEQTTHDETVPSVDDIEGRMDEPKSCDK THTCPPCPAPELLGGPSVFLFPPKPKDTLMISRTPEVTCVWVDVSHEDPEVKFNWYVDGVEVHNAKTKPREE QYNSTYRVVSVLTVLHQDWLNGKEYKCKVSNKALPAPIEKTISKAKGQPREPQVYTLPPSREEMTKNQVSLTC LVKGFYPSDIAVEWESNGQPENNYKTTTPVLDSDGSFFLYSKLTVDKSRWQQGNVFCSSVMHEALHNHYTQ KSLSLSPGK
Formulation	Lyophilized from a 0.2 μm filtered solution of PBS, pH 7.4 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only