

Recombinant Mouse C-X-C Motif Chemokine 16,CXCL16,SR-PSOX (C-6His)

Catalog No: #AP75549

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Package Size: #AP75549-1 10ug #AP75549-2 50ug #AP75549-3 500ug #AP75549-4 1mg

Description

Product Name	Recombinant Mouse C-X-C Motif Chemokine 16,CXCL16,SR-PSOX (C-6His)
Host Species	Human Cells
Purification	Greater than 95% as determined by reducing SDS-PAGE.
Specificity	Mouse
Accession No.	Q8BSU2
Uniprot	Q8BSU2
GeneID	66102;
Calculated MW	20.1kDa
Target Sequence	NQGSVAGSCSCDRTISSGTQIPQGTLDHIRKYLKAFHRCPPFFIRFQLQSKSVCGGSQDQWVRELVDGFERKE CGTGHGKSFHHQKHLQPASTQTPEAAEGTPSDTSTPAHSQSTQHSTLPSGALSLNKEHTQPWEMTTLPSGY GLEARPEAEANEKQDDRQEQEAPGAGASTPAWVDHHHHHH
Formulation	Lyophilized from a 0.2 µm filtered solution of PBS, pH 7.4 .
Storage	Lyophilized protein should be stored at < -20C, though stable at room temperature for 3 weeks.Reconstituted protein solution can be stored at 4-7C for 2-7 days.Aliquots of reconstituted samples are stable at < -20C for 3 months.

Note: This product is for in vitro research use only